


# Spectronon and SpectrononPro Software

Data visualization, analysis, and acquisition

Resonon's  
Window's based  
data visualization  
and acquisition  
software



## Spectronon

Spectronon software provides intuitive visualization and analysis tools for hyperspectral data. Spectronon can be downloaded for free at:  
<http://downloads.resonon.com/>

## SpectrononPro

SpectrononPro provides all the tools of Spectronon, and also enables data acquisition with Resonon imaging spectrometers. SpectrononPro integrates seamlessly with Resonon scanning stages for plug-and-play operation. SpectrononPro is provided with the purchase of all Resonon imaging spectrometers.

## Features

Spectronon provides numerous intuitive tools that enable collecting spectra, classifying images, binning, and much more. Socket connections enable interfacing to Matlab, LabView, or custom built software. Spectronon is fully ENVI® compatible and supports BIL, BIP, and BSQ formatting. For a detailed description, see the Spectronon manual available at:  
<http://www.resonon.com/downloads.html>

## Hardware

For optimal performance, see the hardware recommendations provided at:  
<http://www.resonon.com/spectronon-requirements.html>.